


# MDP2060-0900 for WELCH ALLYN PATIENT MONITORS Large Screen Remote

Displays up to seven waveforms, numerics and alarm information



## Vital Signs Displayable:

- Heart Rate + 2 leads of ECG
- Non-Invasive Blood Pressure
- Invasive Blood Pressure
- Pulse Oximetry
- Impedance Respiration
- Temperature
- Capnography

## Applications:

- Emergency Rooms
- Gastro-Intestinal Endoscopy Departments
- Intensive and Critical Care Units
- Post-Anesthesia Care Units

## Features:

- HIPAA compliant
- Provides visual and audible alarms
- Displays up to seven waveforms

Call VIDCO @ 800-638-4326 to discuss your monitoring requirements.

# MDP2060-0900 Specifications

## MONITOR COMPATIBILITY

Model/Software Rev: Propaq 100 / 8.40.00  
 Propaq Encore / 1.31.01 to 2.40.00  
 Propaq CS / 3.00.00 to 3.50.00

## TIME BASE

25mm/sec & 50mm/sec ( CO2 and RESP which are at 6.25 mm/sec).

## VITAL SIGNS

Waveforms: ECG1, ECG2, RESP, INVP1, INVP2, SPO2, CO2  
 Numeric Values: Heart or Pressure Rate, Blood Pressures, Systolic, Diastolic and Mean values, NIBP, SpO2, EtCO2/InCO2, Resp/Breath Rate, Temperature

## ALARMS

Visual: Flashing yellow bell; Alarming numerics flash red  
 Audible Suspend: 3 minutes when suspended at central station

## TIMER

Start / Stop/Reset

## LANGUAGES

English/German/Spanish/French

## CONTROL

Front Panel Knob or Hand Held Remote

## VIDEO OUTPUT

Analog RGB 0.7Vpp, 75 ohms  
 Sync: TTL Negative going  
Resolution: Horizontal Rate: Vertical Rate:  
 1024 x 768 48.5 kHz 60.16Hz  
 1280 x 1024 64.8 kHz 60.8

## PHYSICAL

	<u>MDP2060-0900/0901</u>	<u>Power Adapter</u>
Size:	H.2.3" (58 mm) W 10.2" (259 mm) D 8.9" (226 mm)	H 1.3" (32 mm) W 2.3" (58 mm) D 4.3" (110 mm)
Weight:	3.0 lbs. (1.36 kg)	0.4 lbs. (0.18 kg)

## AC INPUT

100-240V, 50-60Hz, .4A at 100V, .25A at 240V

## ELECTRICAL SAFETY

Conforms to IEC 601-1-2 as a Class IIa Device  
 Listed for UL2601, CAN/CSA C22.2 NO.601.1

## ENVIRONMENTAL

Temperature  
 Operating: 0 to 40° C  
 Storage: -40 to +85° C  
 Humidity: 5% to 95%, non-condensing  
 EMC: EN55011:1998 (Amended by A1: 1999 and A2: 2002) Class B CISPR 11: 2004 Amended by A2: 2006

## FDA

Device 510K listed

## CE/MEDICAL DEVICE DIRECTIVE

CB Scheme testing, CE marking per Medical Device Directive, Class IIa Device

## QUALITY CERTIFICATION

Certified to the ISO 9001 and ISO 13485 Quality Standards

Vital signs available for display from each Welch Allyn patient monitor

	Quicksigns	Propaq 100 Series	Propaq CS & Encore Series
<b>Waveforms</b>			
ECG1	--	X	X
ECG2	--	--	X
P1	--	X	X
P2	--	X	X
SPO2	--	X	X
CO2	--	X	X
RESP	--	--	X
<b>Numeric Values</b>			
Heart Rate	X	X	X
P1	--	X	X
P2	--	X	X
NIBP	X	X	X
SPO2	X	X	X
CO2	--	X	X
RR (Respiration Rate)	--	BR only	X
BR (Breath Rate from CO2)	--		
TEMP	T1	T1	T1,T2

Propaq, Quicksigns and Atlas are registered trademarks of Welch Allyn, Inc.

